DNA analysis workshop

#1 Extract DNA

A) With a DNA-extraction kit

- recommended for samples with strong cell walls fungi, most plants, seeds, cyanobacteria, etc.
- follow the protocol provided with the kit (this is specific for the organism you want to use)
- this will usually yield much better quality & quantity of DNA than necessary

B) With a DIY protocol

• for soft samples - mammalian cells, many fruits, ground up plant tissue, etc.

#1 mush up the sample in a bag or other container (¼ banana, 3-5 berries, etc.) #2 solve a small spoon of salt into ½ cup of hot water → add 3-5 big spoons to sample → mix for 1 min #3 add a small spoonful of dish-washing soap → mix for 1 min

→ you can try the mush directly for the PCR (it's very dirty but often it works) #4 run the mush through coffee-filter

→ you can try the filtrate directly for the PCR (it's not really pure but often it works) #5 slowly pour 2-3 big spoonful of ice-cold isopropyl-alcohol/isopropanol down the side of the container #6 wait ~8 min and observe the DNA forming cloudy clumps → remove some with a big pipette into 'eppi' #7 centrifuge for 20 sec, remove liquid, wash: add 1 ml water, centrifuge, remove water, mix in 50µl water

 \rightarrow use ~2 μl of this solution for the PCR

C) Without any protocol

- some samples can be smeared directly into the PCR-tube (just try and you will see...)
- use veeeery little, so little that you can almost not see anything on the tip of the pipette
- the PCR program should then start with 10 min at 95 $^{\circ}$ C to boil and break the sample

#2 Prepare Overlord-Mastermix

- this should be enough for 9 samples with 10 μ l each (technically 10 but calculate with a bit extra)
- divide the Overlord-Mastermix into small PCR-eppis and add your samples to each of them
- number the tubes and document in your labbook/notes

| Master-Mix (5x) | 20 µl |
|---------------------|--------|
| Primer #1 | 10 µl |
| Primer #1 | 10 µl |
| Loading buffer (6x) | 17 µl |
| Water | 43 µl |
| End-Volume | 100 µl |

#3 Run the PCR

• place the PCR-eppis into the PCR-machine and run the following program

| | Temp | Time | # of cycles |
|----------------------|------|--------|-------------|
| Initial Denaturation | 94°C | 10 min | |
| Denaturation | 94°C | 30 sec | |
| Primer Annealing | 55°C | 30 sec | 35 |
| Extension | 72°C | 30 sec | |
| Final Extension | 72°C | 5 min | |

#4 Gel-Electrophoresis

Prepare the Agarose-Gel

- mix 0.5 g Agarose in 50 ml 1 x TBE buffer (dilute from 10x)
- heat in microwave 30 sec full, mix, heat again in short intervals until the solution is clear
- let it cool down until you can touch it with your hand
- add 2.5 µl DNA stain green
- pour into gel-chamber (For BentoLab: make sure the rubbers are in place)
- insert comb for 12 wells or more

Prepare TBE running buffer

- prepare 100 ml 1 x TBE buffer (dilute from 10x)
- add 2.5 µl DNA stain green
- remove the rubbers (only for BentoLab) + comb from the gel-chamber
- use the TBE running buffer to cover the gel (until "max"-marking)
- 100 ml will be enough for 2-3 gels

Load the ladder & samples

- the gel should be completely covered in buffer
- gently pipette 10 μl of the DNA-ladder (premixed with loading buffer and dye) into the first and last well of the gel
- gently pipette 10 μ l of the PCR-sample (all of it) into the other wells (document positions)
- if there is some solid or slimy something in the sample, don't pipette it in

Run the gel

- close the box
- run at 50 V for 40 min
- after observation you might want to run again for 10-20 minutes

Observe the result

- remove the lid
- gently place the gel-chamber onto the illuminator
- place the observation-box on top
- switch on the light and look through the hole
- take a picture